



# Panorama: the next generation of NIPT

Non-invasive prenatal screen

Bangkok  
**Cytogenetics**  
Center

65/18 Soi Vibhavadi-Rangsit 16/6 Vibhavadi-Rangsit Road.  
Kwang Chomphol, ChatuchakBangkok 10900 Thailand  
Tel : +662 6900063, Fax : +662 6900064  
[www.bccgroup-thailand.com/panorama](http://www.bccgroup-thailand.com/panorama)

 : Panorama NIPT Thailand



**Panorama™**  
Next-generation NIPT  
POWERED BY NATERA


# Panorama improves upon first-generation NIPTs

## First-generation NIPTs use a “counting methodology” to assess risk

Laboratories that utilize whole-genome sequencing technologies (referred in this brochure as WGS-1 and WGS-2) and array-NIPT, examine fragments of conserved DNA sequences – the 99% of our DNA that makes us the same. These labs compare counts of fragments from chromosomes of interest, such as chromosome 21, against a selected reference chromosome, such as chromosome 3.

If the ratio of fragments between the chromosome of interest and the reference chromosome is determined by the lab to be out of proportion, then the lab identifies the result as “high risk”

### Counting methodology



By looking at conserved DNA sequences and not distinguishing between maternal and fetal DNA, counting methodologies cannot detect triploidy, vanished twin, maternal mosaicism, and complete molar pregnancies.

Failure to identify these conditions can result in false negatives, false positives, and delayed diagnosis of conditions associated with maternal complications.

## Panorama’s SNP-based technology offers greater accuracy than first-generation NIPTs<sup>1-11</sup>

Panorama provides results with fewer false negatives, fewer false positives, and identification of maternal complications.


\*Representation of counting methodology for illustrative purposes

## Panorama is the only NIPT that can distinguish between maternal and fetal (placental) DNA

Panorama isolates single nucleotide polymorphisms (SNPs) – the 1% of our DNA that makes us different from one another.


Our technology sequences targeted chromosomal regions of interest and plots SNP patterns from maternal and fetal cell-free DNA. The patterns are evaluated by our proprietary algorithm to determine if the allele patterns indicate increased risk of fetal abnormalities.

### Panorama's SNP-based methodology



By distinguishing between maternal and fetal DNA, Panorama can detect triploidy, vanished twin, and complete molar pregnancies. This distinction also minimizes the chance that maternal mosaicism will lead to an incorrect result.

### Compared to first-generation NIPTs, Panorama reduces both false negative rates (FNR) and false positive rates (FPR)<sup>12-16</sup>



Array-NIPT is excluded from the FPR chart because data on monosomy X, a leading contributor to false positives for counting methodologies, is not reported in reviewed literature.<sup>4</sup>


\*Representation of a SNP plot for illustrative purposes

## Panorama's SNP-based technology results in the highest validated fetal sex accuracy of any NIPT<sup>1-4,8,9,17</sup>

Panorama utilizes a specific sex-chromosome algorithm that compares SNPs from X and Y to determine the presence and copy number of Y.<sup>18</sup>

With first generation NIPTs, as many as 1 in 77 cases may report incorrect gender. A wrong call can lead to unnecessary clinical work-up and create anxiety for the patient.


### Fetal sex error rates: summary of validation studies



## Panorama's SNP-based approach yields the highest commercially available sensitivity for 22q<sup>19-21</sup>

By evaluating unique DNA sequences within the critical region associated with 22 q11.2 deletion syndrome, Panorama has a higher detection rate than counting methodologies. First generation NIPTs count conserved DNA fragments for chromosome 22 and can overlook small deletions, like 22q.

### Panorama leads the field in 22q11.2 screening sensitivity




## Accurate fetal fraction measurement is essential to accurate results<sup>21</sup>

Panorama is the only NIPT that has always measured and reported fetal fraction.

## Panorama's SNP-based method is a gold standard in fetal fraction measurement


	Panorama <sup>1,2,3</sup>	Array-NIPT <sup>1,2,3</sup>	WGS-1 NIPT <sup>5,6,7,24</sup>	WGS-2 NIPT <sup>9,10,11</sup>
Method of fetal fraction measurement	13,392 SNPs	576 SNPs	Distribution of short (<150 bp) cfDNA	No data available on methodology or performance
Combined false negative rate in validation studies (Trisomies 21, 18, 13)	0.60%	1.33%	1.89%	2.40%

## Counting methodologies' ability to detect abnormalities drops off below 8% fetal fraction, which can produce false negative results<sup>25,26</sup>



## Deeper sequencing on chromosomal regions of interest enables Panorama to maintain high-quality results at lower fetal fractions

Panorama's proprietary algorithm incorporates fetal fraction measurement and reflexes samples with lower fetal fraction to a higher depth of read.



# Are you offering Panorama to women of all ages?

## NIPT is strongly supported by guidelines

The American Congress of Obstetricians and Gynecologists (ACOG), as well as the American College of Medical Genetics and Genomics (ACMG), among other societies, now acknowledges the use of NIPT for all singleton pregnancies, regardless of age or risk.<sup>29,30</sup>

## Panorama is the only NIPT validated in high- and low-risk patients

Validation T21, T18, T13, and MX <sup>2</sup>		
<b>High-risk:</b>	Sensitivity: 98.0%	Specificity: 99.5%
<b>Low-risk:</b>	Sensitivity: 100%	Specificity: 100%

## Professional societies recognize NIPT as a first-line screening option



“Informing all pregnant women that NIPS is the most sensitive screening option”

**ACMG Position Statement, July 2016**



“Data on the performance of cell-free DNA testing in the general obstetric population are now available [and]... similar to the levels previously published for the high-risk population.”

**ACOG/SMFM Practice Bulletin #163, May 2016**



“Different scenarios... are possible, including NIPT as an alternative first-tier option.

**ASHG policy, March 2015**



“The following protocol [is] currently considered appropriate; cfDNA screening as a primary test offered to all pregnant women.”


**Position statement from the Chromosome Abnormality Screening Committee, June 2015**

## Panorama helps clinicians triage twin pregnancies effectively<sup>2-7</sup>

While chorionicity can be reliably detected early in a pregnancy, studies have shown that up to 19% of monozygotic pregnancies are incorrectly classified as dichorionic.<sup>4</sup>

Panorama allows clinicians to align their ultrasound findings with an early and accurate zygosity determination.

Identifying a monozygotic twin pregnancy with Panorama can prompt earlier, targeted ultrasound assessments for chorionicity and associated complications. Knowing that a twin pregnancy is dizygotic reduces concerns about TTTS.



## References

1. Natera validation data. Manuscript in preparation.
2. Society for Maternal-Fetal Medicine, Clinical guideline: Twin-twin transfusion syndrome, Jan 2013.
3. American College of Obstetricians and Gynecologists and Society for Maternal-Fetal Medicine, practice bulletin #169, Oct 2016.
4. Blumenfeld et al. J Ultrasound Med. 2014 Dec;33(12):2187-92.
5. Oldenburg et al. Ultrasound Obstet Gynecol. 2012 Jan;39(1); 39: 69-74.
6. Chasen, Chervenak (2017). Twin pregnancy: Prenatal issues. In T. Post (Ed.), UpToDate. Waltham, Mass.: UpToDate. Retrieved from [www.uptodate.com](http://www.uptodate.com)
7. Cunningham et al. Williams Obstetrics. 24th edition. New York: McGraw-Hill Education, 2014.
8. Boyle et al. BJOG. 2014 Feb; 121(7):809-20. doi: 10.1111/1471-0528.12574.
9. Sparks et al. Obstet Gynecol. 2016 Nov;128(5):1127-33.
10. Rodis et al. Obstet Gynecol. 1990 Dec;76(6):1037-41.
11. <http://ghr.nlm.nih.gov/condition/trisomy-18>.
12. <http://ghr.nlm.nih.gov/condition/trisomy-13>.

**Patient Information**

Patient Name: Jane Doe  
 Date of Birth: 11/08/1975  
 Maternal Age at EDD: 37  
 Gestational Age: 11 weeks/0 days  
 Maternal Weight: N/A  
 Patient ID: P99457  
 Medical Record #: M84555  
 Collection Kit: 254233-2-N  
 Accessioning ID: C47695  
 CaseFile ID: 159466

**Test Information**

Ordering Physician: Dr. Matthew Goodbirth,  
 M.D. (G123456)  
 Clinic Information: Natera, Inc.  
 Additional Reports: N/A  
 Report Date: 02/01/2013  
 Samples Collected: 01/31/2013  
 Samples Received: 02/01/2013  
 Mother Blood



ABOUT THIS SCREEN: Panorama™ is a screening test, not diagnostic. It evaluates genetic information in the maternal blood, which is a mixture of maternal and placental DNA, to determine the chance for specific chromosome abnormalities. The test does NOT tell with certainty if a fetus is affected, and only tests for the conditions ordered by the healthcare provider. A low risk result does not guarantee an unaffected fetus.

**FINAL RESULTS SUMMARY : TWINS**

Result

**HIGH RISK for Trisomy 21**

Zygosity

**Dizygotic****FRATERNAL TWINS**

Fetal Sex

**Male****Female**

Fetal Fraction(s)

**8.3%, 8.4%**

This is a screening test only. Genetic counseling and diagnostic testing for both fetuses should be offered to further evaluate these findings.

Panorama analyzes DNA from the placenta. In some cases placental DNA can differ from that of the fetus. Therefore, even with high risk results, the fetus may be unaffected

**RESULT DETAILS: ANEUP LOIDIES**

Condition tested <sup>1</sup>	Result	Risk Before Test <sup>2</sup>	Risk After Test <sup>3</sup>
Trisomy 21	High Risk	1/152	7/10
Trisomy 18	Low Risk	1/111	<1/10,000
Trisomy 13	Low Risk	1/357	<1/10,000

1. Reporting for Monosomy X, Triploidy, and microdeletion syndromes is not available for dizygotic twin pregnancies. Excludes cases with evidence of fetal and/or placental mosaicism. 2. Based on maternal age, gestational age, and/or general population, as applicable. References available upon request. 3. Risk after test for aneuploidy incorporates results from the Panorama algorithm as well as analytical PPV (high risk) and NPV (low risk). Maternal age is utilized in this calculation, however the "risk after test" may not reflect the actual PPV for this patient, as additional risk factors, including but not limited to: results of other screening, ultrasound findings, personal/family history, are not included in the risk assessment.

Approved By:  Susan Zneimer, Ph.D., FACMG, Laboratory Director

IF THE ORDERING PROVIDER HAS QUESTIONS OR WISHES TO DISCUSS THE RESULTS, PLEASE CONTACT US AT 650-249-9090 #3. Ask for the NIPT genetic counselor on call.



**Patient Information**

Patient Name: Jane Doe  
 Date of Birth: 11/08/1975  
 Maternal Age at EDD: 37  
 Gestational Age: 11 weeks/0 days  
 Maternal Weight: N/A  
 Patient ID: P99457  
 Medical Record #: M84555  
 Collection Kit: 254233-2-N  
 Accessioning ID: C47695  
 Case File ID: 159466

**Test Information**

Ordering Physician: Dr. Matthew Goodbirth,  
 M.D. (G123456)  
 Clinic Information: Natera, Inc.  
 Additional Reports: N/A  
 Report Date: 02/01/2013  
 Samples Collected: 01/31/2013  
 Samples Received: 02/01/2013  
 Mother Blood

**ABOUT THIS SCREEN:** Panorama™ is a screening test, not diagnostic. It evaluates genetic information in the maternal blood, which is a mixture of maternal and placental DNA, to determine the chance for specific chromosome abnormalities. The test does NOT tell with certainty if a fetus is affected, and only tests for the conditions ordered by the healthcare provider. A low risk result does not guarantee an unaffected fetus.

**FINAL RESULTS SUMMARY**

Result

**HIGH RISK for Angelman syndrome**

Fetal Sex

**Male**

Fetal Fraction

**8.3%**

This is a screening test only. Genetic counseling and diagnostic testing, either by microarray and UPD testing for both iso- and heterodisomy, or by methylation testing, should be offered to further evaluate these findings.

Panorama analyzes DNA from the placenta. In some cases placental DNA can differ from that of the fetus. Therefore, even with high risk results, the fetus may be unaffected.

**RESULT DETAILS: ANEUPLOIDIES**

Condition tested <sup>1</sup>	Result	Risk Before Test <sup>2</sup>	Risk After Test <sup>3</sup>
Trisomy 21	Low Risk	<b>1/152</b>	<b>&lt;1/10,000</b>
Trisomy 18	Low Risk	1/354	<1/10,000
Trisomy 13	Low Risk	1/1,116	<1/10,000
Monosomy X	Low Risk	1/255	<1/10,000
Triploidy	Low Risk		

**RESULT DETAILS: MICRODELETIONS**

Condition tested <sup>1</sup>	Result	Risk Before Test <sup>2</sup>	Risk After Test <sup>4</sup>
22q11.2 deletion syndrome	Low Risk	<b>1/2,000</b>	<b>1/9,000</b>
1p36 deletion syndrome	Low Risk	1/5,000	1/12,400
Angelman syndrome	High Risk	1/12,000	1/10
Cri-du-chat syndrome	Low Risk	1/20,000	1/57,100
Prader-Willi syndrome	Risk Unchanged	1/10,000	1/10,000

1. Excludes cases with evidence of fetal and/or placental mosaicism. 2. Based on maternal age, gestational age, and/or general population, as applicable. References available upon request. 3. Risk after test for aneuploidy incorporates results from the Panorama algorithm and data from a published study of 17,885 women [Dar et al. Am J Obstet Gynecol. 2014. Nov;211(5):527.e1-27.e17] and are reported as PPV (high risk) and NPV (low risk). Maternal age is utilized in this calculation, however the "risk after test" may not reflect the actual PPV for this patient, as additional risk factors, including but not limited to; results of other screening, ultrasound findings, personal/family history, are not included in the risk assessment. Risk after test for microdeletion(s) incorporates results from the Panorama algorithm and data from published studies [Martin et al. Clin Genetics. 2017 Jul 11; Wapner R J et al. Am J Obstet Gynecol. 2015 Mar;212 (3):332 .e1-9] and are reported as PPV (high risk) and NPV (low risk). Risk for microdeletions is independent of maternal age. Fetal fraction (FF) is utilized in this calculation. Depending upon FF, in some cases only the paternal allele is evaluated (see page 2). The "risk after test" may not reflect the actual PPV for this patient, as additional risk factors, including but not limited to; results of other screening, ultrasound findings, personal/family history, are not included in the risk assessment.

Approved By:  Susan Zneimer, Ph.D., FACMGG, Laboratory Director

IF THE ORDERING PROVIDER HAS QUESTIONS OR WISHES TO DISCUSS THE RESULTS, PLEASE CONTACT US AT 650-249-9090 #3. Ask for the NIPT genetic counselor on call.



ABOUT THIS SCREEN: Panorama™ is a screening test, not diagnostic. It evaluates genetic information in the maternal blood, which is a mixture of maternal and placental DNA, to determine the chance for specific chromosome abnormalities. The test does NOT tell with certainty if a fetus is affected, and only tests for the conditions ordered by the healthcare provider. A low risk result does not guarantee an unaffected fetus.

Under collaboration between faculty of Medicine Siriraj Hospital, Mahidol University and Bangkok Cytogenetics Center Company Limited.

#### Patient Information

Patient Name:  
Date of Birth:  
Maternal Age at EDD: 39  
Gestational Age: 13 weeks/2 days  
Maternal Weight: 38.6 Kgs  
Patient ID:  
Medical Record #: N/A  
Collection Kit: 5979155-2-N  
Reference ID: N/A  
Accession ID:  
Case File ID:

#### Test Information

Ordering Physician:  
Hospital/Clinic:  
Additional Reports: N/A  
Report Date: 13 Nov 2018  
Samples Collected: 05 Nov 2018  
Samples Received: 07 Nov 2018  
Mother Blood

#### FINAL RESULTS SUMMARY

Result

LOW RISK



Fetal Sex

Female



Fetal Fraction

23.55%



#### RESULT DETAILS: ANEUPLOIDIES

Condition tested <sup>1</sup>	Result	Risk Before Test <sup>2</sup>	Risk After Test <sup>3</sup>
Trisomy 21	Low Risk	1/125	<1/10,000
Trisomy 18	Low Risk	1/330	<1/10,000
Trisomy 13	Low Risk	1/1029	<1/10,000
Monosomy X	Low Risk	1/568	<1/10,000
Triploidy	Low Risk		

1. Excludes cases with evidence of fetal and/or placental mosaicism. 2. Based on maternal age, gestational age, and/or general population, as applicable. References available upon request. 3. Risk after test for aneuploidy incorporates results from the Panorama algorithm and data from a published study of 17,885 women [Dar et al. Am J Obstet Gynecol. 2014. Nov;211(5):527.e1-27.e17] and are reported as PPV (high risk) and NPV (low risk). Maternal age is utilized in this calculation, however the "risk after test" may not reflect the actual PPV for this patient, as additional risk factors, including but not limited to: results of other screening, ultrasound findings, personal/family history, are not included in the risk assessment.

Analyzed By: Rachawalan

Rachawalan Suriyasaengsri, B.Sc.  
Scientist

Approved By: N. Pongvarin

Naravat Pongvarin, M.D., Ph.D.  
Head of Clinical Molecular Pathology Laboratory

IF THE ORDERING PROVIDER HAS QUESTIONS OR WISHES TO DISCUSS THE RESULTS, PLEASE CONTACT US AT 02-690-0063, 086-306-2084 OR SEND US AN INQUIRY TO INFO@BCCGROUP-THAILAND.COM



Case File ID: 181107\_0005  
Bangkok Cytogenetics Center Co., Ltd.  
65/18 Soi Vibhavadi-Rangsit 16/6, Vibhavadi-Rangsit Road, Kwang Chomphol, Khet Chatuchak, Bangkok 10900 Thailand  
Tel: 02-690-0063, 086-306-2084, www.bccgroup-thailand.com, Email: info@bccgroup-thailand.com



Test Panel	Panorama™ Prenatal Basic Panel (Thailand / USA)	Panorama™ Plus 22q11.2 Panel (Thailand / USA)	Panorama™ Plus Microdeletions Panel (Thailand / USA)
Method	Next generation Sequencing using SNPs Technology		
Facility Location	Thailand: Siriraj Hospital / USA: Natera, California		
Abnormality Detection	1. Aneuploidy of chromosome 13, 18, 21, X, Y 2. Triploidy	1. Aneuploidy of chromosome 13, 18, 21, X, Y* 2. Triploidy 3. 22q11.2 deletion	1. Aneuploidy of chromosome 13, 18, 21, X, Y* 2. Triploidy 3. 5 microdeletions - 22q11.2 deletion - Prader Willi Syndrome - Angelman Syndrome - Cri Du Chat Syndrome - 1p36 deletion
Minimum Gestation Age	9 weeks (Recommendation: 12 weeks)		
Sample Requirement	Maternal blood 16-20 ml in Streck tubes		
Turn Around Time	10 -14 working days		
Confirmation Test	Karyotype + QF-PCR ( up to 20,000 baht )	Karyotype + Array CGH ( up to 20,000 baht )	
Guarantee (False Negative)	3,500,000 THB (Trisomy 13, 18, 21)		
Twin pregnancy	USA only	USA only (Monozygotic twins)	N / A
Egg donor/ Surrogate on singleton pregnancy	USA only	N / A	N / A

\* Twin pregnancies report sex chromosome aneuploidies only for monozygotic twins

POWERED BY NATERA



## References

1. Nicolaides et al. Prenat Diagn. 2013 June;33(6):575-9.
2. Pergament et al. Obstet Gynecol. 2014 Aug;124(2 Pt 1):210-8.
3. Ryan et al. Fetal Diagn Ther. 2016;40(3):219-223.
4. Stokowski et al. Prenat Diagn. 2015 Oct; DOI: 10.1002/pd.4686.
5. Palomaki et al. Genet Med. 2011 Nov;13(11):913-20.
6. Palomaki et al. Genet Med. 2012 Mar;14(3):296-305.
7. Porreco et al. Am J Obstet Gynecol 2014;210.
8. Mazloom et al. Prenat Diagn 2013;33:591-7.
9. Sehnert et al. MolecularDiagn and Gene 2011.
10. Bianchi et al. Obstet Gynecol. 2012 May;119(5):890-901.
11. Bianchi et al. N Engl J Med 2014;370:799-808.
12. Nicolaides et al. Fetal Diagn Ther. 2014;35(3):212-7.
13. Curnow et al. Am J Obstet Gynecol. 2015 Jan;212(1):79.e1-9.
14. Futch et al. Prenat Diagn 2013;33:569-74.
15. Simon et al. Ultrasound Obstet Gynecol 2015; 46(4):506-7.
16. Wang et al. Clinical Chemistry 60:1, 251–259, 2014.
17. Verinata white paper. Analytical validation of the Verifi prenatal test. 2012
18. Samango-Sprouse et al. Prenat Diagn. 2013;33:1–7.
19. Wapner et al. Am J Obstet Gynecol. 2014; DOI: 10.1016/j.ajog.2014.11.041.
20. Hegelson et al. Prenatal Diagnosis. 2015, 35, 1–6.
21. Commercial protocol not validated; Illumina marketing materials cite "Srinivasan et al. Am J Hum Genet. 2013 Feb 7; 92(2): 167–176" which does not match number of reads used in commercial testing.
22. American College of Obstetricians and Gynecologists (ACOG/SMFM), #640, Sept 2015.
23. Juneau et al. Fetal Diagn Ther. 2014;36(4):282-6.
24. Kim et al. Prenatal Diagnosis 2015, 35, 810–815.
25. Canick, et al. Prenatal Diagnosis 2013, 33, 1 –8.
26. Wright et al. Ultrasound Obstet Gynecol 2015; 45: 48–54.
27. Verifi marketing materials, 2016.
28. Internal data, Natera
29. American College of Obstetricians and Gynecologists (ACOG/SMFM), #163, May 2016.
30. American College of Medical Genetics and Genomics (ACMG), Position Statement, Jul 2016.
31. Palomaki, et al. Genetics in Medicine 2017; DOI:10.1038/gim.2016.194.
32. K Dahl, et al. Ultrasound Obstet Gynecol 2011;38:145.

## CAP accredited and CLIA certified, ISO 13485.

CAUTION: The law restricts these tests/products to sale by or on the order of a physician. Indications, contraindications, warning and instructions for use can be found in the Test/Product labeling supplied with each test/product. Information for use only in countries with applicable health authority test/product registrations.

CAUTION: The Panorama™ Test may be sold only by or on the order of a physician. Not for distribution in the United States.

©2017 Natera, Inc. All Rights Reserved. Natera International, Inc., San Carlos, California, USA. Panorama is a trademark of Natera, Inc. All cited trademarks are the property of their respective owners. Printed in USA. PAN\_MD\_BR\_17\_06\_30\_NAT801513\_INTL



201 Industrial Road, Suite 410 | San Carlos, CA 94070  
1-650-249-9090 | Fax 1-650-730-2272 | [www.natera.com](http://www.natera.com)